# measComp

unknown

2023-July-31

# CONTENTS

1	Required Modules 3											
2	Table of Contents       S         2.1       Overview       S											
	2.2	2.2 Driver for Multi-Function Devices										
		2.2.1	Introduction	6								
		2.2.2	Supported models	7								
		2.2.3	Configuration	29								
		2.2.4	Databases	30								
		2.2.5	Box for USB-CTR08, USB-3104, and USB-1808X	39								
		2.2.6	Box for USB-2408-2AO	41								
		2.2.7	Performance measurements	41								
	2.3	Driver	for the USB-CTR08	48								
		2.3.1	Introduction	48								
		2.3.2	Configuration	53								
		2.3.3	Databases	55								
		2.3.4	Wiring to BCDA BC-020 LEMO Breakout Panels	59								
		2.3.5	Performance measurements	65								
		2.3.6	Restrictions	65								
				-								

# author

Mark Rivers, University of Chicago

This package provides EPICS drivers for the some of the USB and Ethernet I/O modules from Measurement Computing. The software is located in the measComp github repository.

# CHAPTER

# ONE

# **REQUIRED MODULES**

Required module	Required for
EPICS base	Base support
asyn	Driver and device support
autosave	Save/restore support
busy	Busy record support
mca	mca record support
scaler	Scaler record support.
seq	State notation language sequencer. Used in MCS mode with USB-CTR08 and for std.

The required versions of each of the above modules for a specific release of measComp can be determined from the measComp/configure/RELEASE file.

# CHAPTER

TWO

# TABLE OF CONTENTS

# 2.1 Overview

This package provides EPICS drivers for the some of the USB and Ethernet I/O modules from Measurement Computing. Currently the USB-CTR04/08, and multi-function modules (E-1608, USB-1208LS, USB-1208FS, USB-1608G, USB-1608GX-2AO, USB-1808/1808X, USB-231, USB-2408-2AO, E-TC, TC-32, USB-TEMP, USB-TEMP-AI, and E-DIO24) are supported. The multi-function modules support analog input and/or output, temperature input (USB-2408-2AO, USB-TEMP, USB-TEMP-AI, E-TC, TC-32), digital input/output, pulse counters (all but TC-32), and pulse generators (USB-1608G and USB-1608GX-2AO).

Support for other modules is straightforward to add and can be done as the demand arises.

This module is supported on both Windows and Linux, 64-bit and 32-bit.

On Windows it uses the Measurement Computing "Universal Library" (UL), which is only available on Windows.

In R4-0 and later it uses the UL for Linux library from Measurement Computing for Linux drivers. This is an [open-source library available on Github](https://github.com/mccdaq/uldaq). The Linux Universal Library API is similar to the Windows UL API, but the functions have different names and different syntax.

UL for Windows and Linux support most current Measurement Computing models.

In versions prior to R4-0 the Linux support used the [low-level drivers from Warren Jasper](https://github.com/wjasper/ Linux\_Drivers). On top of these drivers the module provides a layer that emulates the Windows UL library from Measurement Computing. The EPICS drivers thus always use the Windows UL API and are identical on Linux and Windows. The E-1608, E-TC, E-TC32, E-DIO24, USB-1608G-2AO, USB-CTR08, USB-TEMP, USB-TEMP-AI and USB-31XX models are supported in these versions.

# 2.2 Driver for Multi-Function Devices

#### author

Mark Rivers, University of Chicago

#### Contents

- Driver for Multi-Function Devices
  - Introduction
  - Supported models
    - \* *E-1608*



# 2.2.1 Introduction

This is an EPICS driver for the multi-function devices from MeasurementComputing. These multi-function devices support support analog input, temperature input (thermocouple, RTD, thermistor, and semiconductor), analog output, binary I/O, counters, and timers. Not all devices have all of these capabilities.

The driver is written in C++, and consists of a class that inherits from asynPortDriver, which is part of the EPICS asyn module.

The driver is written to be general, so that it can be used with any Measurement Computing multi-function module. It uses the introspection capabilities of their UL library to query many of the device features. However, there are some features that cannot be queried, so the driver does require small modifications to be be used with a new model.

# 2.2.2 Supported models

The following models are currently supported.

#### E-1608



## Fig. 1: Photo of E-1608

This module costs \$525 and has the following features:

- 16-bit analog inputs
  - 8 single-ended channels or 4 differential channels
  - Programmable per-channel range: +-1V, +-2V, +-5V, +-10V
  - 250 kHz total maximum input rate, i.e. 1 channel at 250 kHz, 2 channels at 125 kHz, etc.
  - Internal or external trigger.
  - Internal or external clock for input signals.
  - Input FIFO, unlimited waveform length
- 16-bit analog outputs

- 2 channels, fixed +-10V range
- No output waveform capability
- Digital inputs/outputs
  - 8 signals, individually programmable as inputs or outputs
- Counter
  - 1 input
  - 10 MHz maximum rate, 32-bit register

More information can be found in the E-1608 product description.

The following is the main medm screen for controlling the E-1608.

# E-TC

This module costs \$505 and has the following features:

- Ethernet interface.
- 8 thermocouple inputs
  - 8 channels with cold-junction compensation. Types J, K, T, E, R, S, B, and N.
  - 4 samples/s.
- Digital inputs/outputs
  - 8 signals, individually programmable as inputs or outputs
- Counters
  - 1 input
  - 10 MHz maximum rate, 32-bit register

More information can be found in the E-TC product description.

The following is the main medm screen for controlling the E-TC.

# TC-32

This module costs \$1999 and has the following features:

- USB and Ethernet interfaces, either can be used.
- 32 thermocouple inputs
  - 32 channels with cold-junction compensation. Types J, K, T, E, R, S, B, and N.
  - 3 samples/s if reading all 32 channels, faster if reading fewer.
- Digital inputs
  - 8 digital inputs, switch-selectable pullup resistor
- Digital outputs
  - 32 digital inputs, switch-selectable pullup resistor



Fig. 2: E1608\_module.adl



Fig. 3: Photo of E-TC



Fig. 4: ETC\_module.adl



Fig. 5: Photo of TC-32

- Each output can either be controlled by software or can be controlled by the alarm status of the corresponding thermocouple. Flexible alarm configuration, i.e. hysteresis.

More information can be found in the TC-32 product description.

The following is the main medm screen for controlling the TC-32.

X TC32_m	X TG2_module_adl@corvete – 🗆 X																				
						TC	-32	TC32	:												
	Model TC32 Unique					)	UL ve	rsion	1.2.0		Poll	sleep	time	(ms) 🆻	0.0						
Model # 306 Firmwar			re 1.08		Dri	ver ve	rsion	4.3		Poll	cycle	time	(ms)5	1.1							
	Thermocouple 0-7		Thermo	couple	e 8-15	5			The	rmocc	uple	16-2	3			The	rmoco	uple	24-31		
0	25.12 1 second = Read	8	29.98 1 se	cond 💷	Read			16	29,6	1 sec	ond II	Read			24	31, 2	1 seco	and 🖃	Read		
1	29,20 1 second a Read	9	29.94 1 se	cond 💷	Read			17	29, 9	9 1 sec	ond 💷	Read			25	29, 14	1 seco	and 💷	Read		
2	30,58 1 second a Read	10	29.48 1 se	cond 💷	Read			18	30.3	3 1 sec	ond 💷	Read			26	30,59	1 seco	and I	Read		
3	20, U3 1 second I Read		29.43 1 se	cond =	Read			19	29, 2	1 sec	ond II	Read		-1	27	31, 19	1 seco	and I	Read		
4	31,12 1 second a Read	12	20.17	cond II	Read			20	20.1	1 1 500	ond II	Read		-1	28	29,20	1 seco	and	Read		
6	30.45 1 second at Read	14	29.87 1 5	cond a	Read			22	30.3	9 1 500	and al	Read		-1	20	30.29	1 5000	and al	Read		5B
7	29,56 1 second at Read	15	28.92 1 54	cond a	Read			23	29.6	9 1 sec	and II	Read		-1	31	29.9	1 seco	and II	Read		- H
	BStrip charts BConfigure		ĐStrip cha	rts 🛛 🖸	Configure	2			PiStr	ip chart	s Pac	onfigure	•			PuStr	ip chart	s Qu	onfigure		
	Digital Inputs										Dig	ital	Outpu	its							_
0	1 2 3 4 6 5	7		0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15		
High	High High High High High High High OXII					Low	Low	Low 0	x0												
		Low High	Low High	Lov High	Low High	Low High	Lov High	Low High	Lou High	Lov High	Low High	Low High	Lov High	Low High	Low High	Low High	Low High 0:	ж0			
		16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31				
				Low	Low	Low															
				Low High	Lou High	Lov High	Low High	Lou High	Lov High	Low High											



#### USB-1608G and USB-1608GX-2AO

This module costs \$799 and has the following features:

- 16-bit analog inputs
  - 16 single-ended channels or 8 differential channels
  - Programmable per-channel range: +-1V, +-2V, +-5V, +-10V
  - 500 kHz total maximum input rate, i.e. 1 channel at 500 kHz, 8 channels at 62.5 kHz, etc.
  - Internal or external trigger. External trigger shared with analog outputs.
  - Internal or external clock, input and output signals.



## Fig. 7: Photo of USB-1608GX-2AO

- 4 kSample input FIFO, unlimited waveform length
- 16-bit analog outputs
  - 2 channels, fixed +-10V range
  - 500 kHz total maximum output rate, i.e. 1 channel at 500 kHz, 2 channels at 250 kHz
  - Internal or external trigger. External trigger shared with analog inputs.
  - Internal or external clock, input and output signals
  - 2 kSample output FIFO, unlimited waveform length
- Digital inputs/outputs
  - 8 signals, individually programmable as inputs or outputs
- Pulse generator
  - 1 output
  - 64MHz clock, 32-bit registers
  - Programmable period, width, number of pulses, polarity
- Counters
  - 2 inputs
  - 20 MHz maximum rate, 32-bit registers

More information can be found in the USB-1608GX-2AO product description.

The USB-1608G is very similar to the USB-1608GX-2AO except that it does not have any analog outputs and the analog inputs are limited to 250 kHz rather than 500 kHz. More information can be found in the USB-1608G product description.

X USB1608G_2AO_module.adl@corvette		– 🗆 X
USB-1608G-2A	0 USB1608G_2A0:	
Model USB-1608GXUnique ID 016669C2UModel #274Firmware 1.09Drive	L version 1.2.0 Poll slee r version 4.3 Poll cycl	p time (ms) <mark>50.0</mark> e time (ms) 58.7
Trigger Mode Positive edge December 2 0 Reset	Digital 1 2 3 4 5 Inputs O O O O	I/0 6 7 8 ● ● 0x0 Low Low 0x0
WaveTorm generator Trigger Internal Current point 0	Direction In In In In In Out Out Out Out Out	High High High In In In Out Out
Retrigger       Disable       # points       1000         Trigger count       Frequency       2.0000         Clock       Internal       Time/point       0.0000         Repeat       Continuous       Total time       0.0000         Output 1       Output 2       Enable       Disable       #         Waveform       Sin wave       Maveform       User-defined       Amplitude #.0000       0ffset 0.0000       0ffset 0.0000       Pulse width 0.0010       Pulse width 0.0010       Pulse width 0.0010       Pulse width 0.0010       Pre-defined waveforms	Analog input           1         0.9999         I/0 Intr         Read           2         5.9478         1 second         Read           3         2.9873         1 second         Read           4         2.9751         1 second         Read           5         1.1957         1 second         Read           6         1.1888         1 second         Read           7         0.5942         1 second         Read           8         0.6888         1 second         Read           Mode         Single-ended         PDonfigure	Waveform digitzer Current point 0 # points 2000 Time/point 0.0000 Total time 0.0000 Time/point 0.0010 # points 2000 First chan 1 # chans 8 Burst mode Disable Trigger Internal
Frequency 2.000 Time/point 0.0010 # points 1 Start Stop Done GPLots	Analog output 1 1.0000 > 2 0.0000 > <u>QConfigure</u>	Retrigger Disable # Trigger count 1 Clock Internal # Continuous One-shot # Auto restart Disable #
Pulse generator           Frequency         Period         Duty Cycle         Width         Delay           5,0000         2,0000e-01         0.5000         1,0000e-01         0.0000e+00           1         5.0000         2.0000e-01         0.5000         1,0000e-01         0.0000e+00	0 State Honorom 1 December 1997	Read rate 1 second # Read Read Start Star Done PPlots

The following is the main medm screen for controlling the USB-1608GX-2AO.

Fig. 8: 1608G\_module.adl

# USB-1808 and USB-1808X

These modules cost \$769 and \$989 and have the following features:

- 18-bit analog inputs
  - 8 single-ended or differential channels
  - Programmable per-channel range:+-5V, +-10V, 0-5V, 0-10V
  - USB-1808: 125 kHz total maximum input rate, i.e. 1 channel at 125 kHz, 8 channels at 15.625 kHz, etc.
  - USB-1808X: 500 kHz total maximum input rate, i.e. 1 channel at 500 kHz, 8 channels at 62.5 kHz, etc.
  - Internal or external trigger. External trigger shared with analog outputs.
  - Internal or external clock, input and output signals.
  - 4 kSample input FIFO, unlimited waveform length
- 16-bit analog outputs
  - 2 channels, fixed +-10V range



Fig. 9: Photo of USB-1808

- USB-1808: 250 kHz total maximum output rate, i.e. 1 channel at 250 kHz, 2 channels at 125 kHz
- USB-1808X: 1000 kHz total maximum output rate, i.e. 1 channel at 1000 kHz, 2 channels at 500 kHz
- Internal or external trigger. External trigger shared with analog inputs.
- Internal or external clock, input and output signals
- 2 kSample output FIFO, unlimited waveform length
- Digital inputs/outputs
  - 4 signals, individually programmable as inputs or outputs
- · Pulse generator
  - 2 outputs
  - 100 MHz clock, 32-bit registers
  - Programmable period, width, number of pulses, polarity
- Counters
  - 2 inputs
  - 50 MHz maximum rate, 32-bit registers
- Quadrature encoder inputs
  - 2 inputs
  - 50 MHz maximum rate, 32-bit registers

More information can be found in the USB-1808 product description.

The following is the main medm screen for controlling the USB-1808.

#### USB-2408-2AO

This module costs \$699 and has the following features:

- 24-bit analog inputs
  - 16 single-ended channels or 8 differential channels
  - Programmable per-channel range: 8 ranges from +-0.078V to +-10V
  - Thermocouple support for 8 channels with cold-junction compensation. Types J, K, T, E, R, S, B, or N.
  - 1 kHz total maximum input rate, i.e. 1 channel at 1 kHz, 8 channels at 125 Hz, etc.
  - Input FIFO, unlimited waveform length
- 16-bit analog outputs
  - 2 channels, fixed +-10V range
  - 1000 Hz total maximum output rate, i.e. 1 channel at 1000 Hz, 2 channels at 500 Hz
  - Output FIFO, unlimited waveform length
- Digital inputs/outputs
  - 8 signals, individually programmable as inputs or outputs
- Counters
  - 2 inputs



Fig. 10: 1808\_module.adl



Fig. 11: Photo of Photo of USB-2408-2AO

- 1 MHz maximum rate, 32-bit registers

More information can be found in the USB-2408-2AO product description. The following is the main medm screen for controlling the USB-2408-2AO.

X USB2408_module.adl@corvette			- 🗆 X
	USB-2408-2A0	USB2408:	
Model USB-2408-2 Unic Model <b>#</b> 254 Fir	ue ID 01AAA83E UL ve mware1.01 Driver ve	ersion 1.2.0 Poll slee ersion 4.2 Poll cycl	p time (ms) <mark>50.0</mark> .e time (ms) 62.8
Analog input 1 0.0000 1 second # Read 2 0.0000 1 second # Read 3 0.0000 1 second # Read 4 0.0000 1 second # Read 5 0.0000 1 second # Read	Counters 1 0 Reset 2 0 Reset C Dir	Digital 1 2 3 4 5 Inputs Low Low Low Low Low High High High High rection In In In In Out Out Out Out Out	I/O 6 7 8 • • • Oxff I.ov Low Low 0xff I.in In In Out Out Out
6 0.0000 1 second # Read 7 0.0000 1 second # Read 8 0.0000 1 second # Read 8 0.0000 1 second # Read 9 EStrip charts @Donfigure 1 27.02 1 second # Read 2 27.30 1 second # Read 3 0.00 1 second # Read 4 0.00 1 second # Read 5 0.00 1 second # Read 6 0.00 1 second # Read 6 0.00 1 second # Read 7 0.00 1 second # Read 8 0.00 1 second # Read 8 0.00 1 second # Read	Waveform Trigger Internal Retrigger Disable Trigger count 1 Clock Internal Repeat One-shot Output 1 Enable Disable Waveform User-defined Amplitude 1.0000 Offset 0.0000 Pulse width 0.0010 User-defined waveforms Frequency 0.0000	generator Current point 0 # points 1000 Frequency 0.0000 Time/point 0.0000 Total time 0.0000 Output 2 Enable Disable Amplitude 1.0000 Offset 0.0000 Pulse width 0.010 Pre-defined waveforms Frequency 0.0000	Waveform digitzer Current point 0 # points 1 Time/point 0.0000 Total time 0.0000 Time/point 0.0010 # points 1 First chan 1 # chans 1 Burst mode Disable Trigger Internal Retrigger Disable Trigger count 1 Clock Internal Continuous One-shot T
Analog output	Time/point inf # points 1000 Start Stop Done	Time/point inf # points 1000 PaPlots	Auto restart Disable Read rate <u>1 second</u> Read Read Start Stop Done CPlots

Fig. 12: 2408\_module.adl

# **USB-TEMP and USB-TEMP-AI**

The USB-TEMP costs \$605 and the USB-TEMP-AI costs \$795. They have the following features:

- Temperature inputs
  - 8 temperature inputs on USB-TEMP, 4 on USB-TEMP-AI. These can be platinum resistance thermometers (RTD), thermocouples, thermistors, or semiconductor sensors.
  - Thermocouple support has cold-junction compensation. Types J, K, T, E, R, S, B, or N.
  - 2 samples/s per channel.
- 24-bit analog inputs (USB-TEMP-AI only)
  - 4 channels
  - Programmable per-channel range: 4 ranges from +-1.25V to +-10V
- Digital inputs/outputs



Fig. 13: Photo of Photo of USB-TEMP

- 8 signals, individually programmable as inputs or outputs
- Counters
  - 1 input
  - 1 MHz maximum rate, 32-bit register

More information can be found in the USB-TEMP product description.

The USB-TEMP and USB-TEMP-AI behave differently from all other Measurement Computing devices. On Windows InstaCal is used to select the temperature sensor type (RTD, thermocouple, etc.) and the RTD wiring configuration. Those settings are written into non-volatile memory on the device, and cannot be changed with EPICS. However, they **can** be changed with EPICS on Linux, so they are exposed in the OPI screen.

The following is the main medm screen for controlling the USB-TEMP-AI.



Fig. 14: USBTEMP\_AI\_module.adl

× measCompUSBTempSetup4.adl@corvette Temperature setup USB\_TEMP\_AI: Display TC Type Scale Sensor Wiring Min Max Ti1 RTD Celsius 🖬 0.00 0.00 2 wire 2 sensors Type J Ti2 RTD 2 wire 2 sensors Celsius **□** 0.00 0,00 Type J **⊒|0.00** 0,00 Ti3 Disabled 2 wire 2 sensors Celsius Type J **= 0.00** 0.00 Ti4 Disabled 2 wire 2 sensors 1 Type J Celsius

The following is the screen for configuring the temperature inputs.



#### **USB-1208LS**

This module costs \$129 and has the following features:

- 12-bit analog inputs
  - 4 differential channels
  - Programmable per-channel range: 8 ranges from +-1V to +-20V
  - 50 Hz maximum sampling rate. The module has a trigger input that allows higher sampling rates, but this
    is not yet supported in the EPICS driver.
- 10-bit analog outputs
  - 2 channels, fixed 0 to +5V range
  - 100 Hz maximum input rate
- Digital inputs/outputs
  - 16 signals, programmable as inputs or outputs in groups of 8
- Counters
  - 1 input
  - 1 MHz maximum rate, 32-bit register

More information can be found in the USB-1208LS product description.

The USB-1208HS, USB-1208FS-Plus and USB-231 are similar devices but with higher performance. These are also supported.

The following is the main medm screen for controlling the USB-1208LS.



Fig. 16: Photo of USB-1208LS



Fig. 17: USB1208LS\_module.adl

# E-DIO24



# Fig. 18: Photo of E-DIO24

This module costs \$320 and has the following features:

- Digital inputs/outputs
  - 24 signals, individually programmable as inputs or outputs
- Counters
  - 1 input
  - 10 MHz maximum rate, 32-bit register

More information can be found in the E-DIO24 product description.

The following is the main medm screen for controlling the E-DIO24.



Fig. 19: EDIO24\_module.adl

## **USB-3100**



#### Fig. 20: Photo of USB-3101

This series of module costs from \$330 (USB-3101) to \$660 (USB-3106) depending on the number of channels and the output type, and has the following features:

- 16-bit analog outputs
  - 4, 8 or 16 channels, individually programmable range 0-10V or +-10V.
  - Some models provide 0-20 mA current output as well as voltage output
  - Some models have high-drive voltage output (+-40 mA)
  - 100 Hz maximum output rate
- Digital inputs/outputs
  - 8 signals, individually programmable as inputs or outputs
- Counters
  - 1 input

- 1 MHz maximum rate, 32-bit register

More information can be found in the USB-3100 series product description.

The following is the main medm screen for controlling the USB-3104 8-channel unit.



Fig. 21: USB3104\_module.adl

The following is the medm screen for configuring the analog outputs on the USB-3104 8-channel unit.

Xn	X measCompAoSetup8.adl@corvette - 🗆 X											
	Analog output setup 13IDC:USB3104:											
	Mode	Range	Drive	limits	Disp	olay	Eng.	units				
1	Normal 🖃	+= 10V 💷	-10,0000	10,0000	-10,0000	10,0000	-10,0000	10.0000				
2	Normal 🗖	+= 10V 💷	-10,0000	10,0000	-10,0000	10,0000	-10,0000	10.0000				
3	Normal 🗖	+= 10V 💷	0.0000	5.0000	-10,0000	10,0000	-10,0000	10.0000				
4	Normal 🗖	+= 10V 💷	0.0000	5.0000	-10,0000	10,0000	-10,0000	10,0000				
5	Normal 🗖	0-10V 🖃	0.0000	10,0000	-10,0000	10.0000	0.0000	10,0000				
6	Normal 🗖	0-10V 🖃	0.0000	10,0000	-10,0000	10.0000	0.0000	10,0000				
7	Normal 🗖	0-10V 🖬	0.0000	5.0000	-10,0000	10.0000	0.0000	10,0000				
8	Normal 🖃	0-10V 🖃	0.0000	5.0000	-10,0000	10,0000	0.0000	10.0000				

Fig. 22: USB3104\_setup.adl

# 2.2.3 Configuration

The following lines are needed in the EPICS startup script for the multifunction driver.

The uniqueID is a string that identifies the device to be controlled.

- For USB devices the uniqueID is the serial number, which is printed on the device (e.g. "01F6335A").
- For Ethernet devices the uniqueID can either be the MAC address (e.g. "00:80:2F:24:53:DE"), or the IP address (e.g. "10.54.160.63", or the IP DNS name (e.g. "gse-e1601-1"). The MAC address, IP address or IP name can be used for devices on the local subnet, while the IP address or IP name must be used for devices on other subnets.

The measComp module comes with example iocBoot/ directories that contain example startup scripts and example substitutions files for each supported model.

# 2.2.4 Databases

The following tables list the database template files that are used with the multi-function modules.

# **Overall Device Functions**

These are the records defined in measCompDevice.template. This database is loaded once for each module.

EPICS	EPICS	asyn in-	drvInfo	Description
record	record	terface	string	
name	type			
\$(P)Mode	INstarinegin	asyn-	MODEL_	NAME model name of this device, e.g. "USB-1808X".
		Octe-		
		tRead		
\$(P)Mode	1 Nangber	asynInt32	MODEL_	NCIM Brokel number of this device, e.g. 318.
\$(P)Firmv	va <b>str¥kegsii</b> on	asyn-	FIRMWA	REht/ERSION version, e.g. "1.03".
		Octe-		
		tRead		
\$(P)Uniqu	ekDringin	asyn-	UNIQUE	IDhe unique ID of this device, e.g. "02151405"
		Octe-		
		tRead		
\$(P)ULVe	rs <b>itmi</b> ngin	asyn-	UL_VERS	STORE version of the UL library on Linux or Windows, e.g. "1.2.0".
		Octe-		
		tRead		
\$(P)Drive	r Vsetrisningnin	asyn-	DRIVER	VEReSkedston of the EPICS driver, e.g. "4.3".
		Octe-		
		tRead		
\$(P)PollT	maiMS	asyn-	POLL_TI	MEheMaStual time for the last poll cycle in ms.
		Float64		
\$(P)PollS	eapoMS	asyn-	POLL_SL	EERe_Mifse to sleep at the end of each poll cycle in ms.
	_	Float64		
\$(P)LastE	rr <b>waMe</b> ssage	asyn-	LAST_ER	RDR_1MESSACOFEssage from the driver.
	form	Octe-		
		tRead		

The medm sub-screen that displays these records. The main screen for every module contains a subscreen like this.

X USB18	808_module.adl@corvette			-	- [	) X
		USB-	1808 13IDC:USB1808:			
	Model USB-1808X	Unique ID 02151405	UL version 1.2.0	Poll sleep time (ms) 50.0		
	Model <b>#</b> <mark>318</mark>	Firmware 1.03	Driver version 4.2	Poll cycle time (ms) 55.7		

Fig. 23: measCompDevice.adl

# Analog I/O Functions

These are the records defined in measCompAnalogIn.template. This database is loaded once for each analog input channel

EPICS	EPICS	asyn in-	drvInfo	Description						
record	record	terface	string							
name	type									
\$(P)\$(R)	ai	asynInt32	ANA-	Analog input value. This is converted from the 16-bit unsigned inte-						
			LOG_IN_	Vedrutevice units from the driver to engineering units using the EGUL						
				and EGUF fields. This value is polled in the driver at the polling fre-						
				quency set by PollSleepMS. The asynInt32Average device support						
				is used, so that the ai value is the average of all the readings from						
				the poller since the last time the record processed. For example,						
				if the poller is running at 100 Hz and the ai record SCAN field is						
				"0.2 seconds" then 20 values will be averaged each time the record						
				processes. If SCAN=I/O Intr then the device support will average						
				the number of values specified in the SVAL field of the record. If						
				SVAL<=1 then the record will processes on each callback, so there						
				is no averaging.						
\$(P)\$(R)R	a <b>ngle</b> bo	asynInt32	ANA-	Input range for this analog input channel. Choices are determined at						
			LOG_IN_	RrANGEne based on the model in use.						
\$(P)\$(R)T	ypnabbo	asynInt32	ANA-	Input type (e.g. "Volts", "TC deg", etc.) for this analog input chan-						
			LOG_IN_	TMAPEChoices are determined at run time based on the model in use.						

The following is the medm screen for controlling the analog input records for the USB-1608GX-2AO. Note that the engineering units limits (EGUL and EGUF) do not have to be in volts, they can be in any units such as "percent", "degrees", etc.

× m	measCompAiSetup.adl												
	Analog input setup USB2408:												
	Mode Range Display Eng. units												
1	TC deg.	-	+= 0,312V 🖬	-10,0000	10.0000	-1,2500	1,2500						
2	Volts	-	+= 0,078V 🖬	-0,1000	0.1000	-0.0780	0.0780						
3	Volts	-	+= 2,5V 🖬	-5,0000	5.0000	-5,0000	5.0000						
4	Volts	-	+= 0,078V 🖬	-5,0000	5.0000	-5,0000	5.0000						
5	Volts	-	+= 5V 🗖	-2,0000	2.0000	-2,0000	2,0000						
6	Volts	-	+= 5V 🗖	-2,0000	2.0000	-2,0000	2,0000						
7	Volts	-	+= 10V 🖃	-1,0000	1.0000	-1,0000	1.0000						
8	Volts	-	+= 10V 💷	-1.0000	1.0000	-1.0000	1,0000						
	-				_								

Fig. 24: measCompAiSetup.adl

These are the records defined in measCompAnalogOut.template. This database is loaded once for each analog output channel

EPICS	EPICS	asyn in-	drvInfo	Description
record	record	terface	string	
name	type			
\$(P)\$(R)	ai	asynInt32	ANA-	Analog output value. This is converted from engineering units to the
			LOG_OU	T10AbitUm signed integer device units for the driver using the EGUL
				and EGUF fields.
\$(P)\$(R)R	la <b>ngke</b> bo	asynInt32	ANA-	Output range for this analog output channel. Choices are determined
			LOG_OU	T <u>a</u> RranGine based on the model in use.
(P)(R)R	le <b>tai</b> rn	asynInt32	ANA-	Analog output value to return to at the end of a pulse. This is con-
			LOG_OU	TvWAddLfFom engineering units to the 16-bit unsigned integer device
				units for the driver using the EGUL and EGUF fields.
\$(P)\$(R)P	u <b>lse</b>	N.A.	N.A.	Choices are "Normal" and "Pulse". In Normal mode the Return
				record is ignored. In Pulse mode the $(P)(R)$ output is written to
				to hardware, followed immediately by writing the \$(P)\$(R)Return
				value.
\$(P)\$(R)T	w <b>aa</b> kVal	N.A.	N.A.	The amount by which to tweak the out when the Tweak record is
				processed.
\$(P)\$(R)T	weakdopt	N.A.	N.A.	Tweaks the output up by TweakVal.
\$(P)\$(R)T	weakDotwn	N.A.	N.A.	Tweaks the output down by TweakVal.

The following is the medm screen for controlling the analog output records for the USB-1608GX-2AO. Note that the engineering units limits (EGUL and EGUF) do not have to be in volts, they can be in any units such as "percent", "degrees", etc. The drive limits can be more restrictive than the full +-10V output range of the analog outputs.

🗙 measComp	AoSetup.adl				
	Analog	outpu	t setup 16	508G:	
Mode	Drive	limits	Display	Eng.	units
1 Normal	-10,0000	10.0000	-10,0000 10,000	0 -10.0000	10,0000
2 Normal	-10,0000	10,0000	-10,0000 10,000	0 -10.0000	10.0000

Fig. 25: measCompAoSetup.adl

# **Temperature Functions**

These are the records defined in measCompTemperatureIn.template. This database is loaded once for each temperature input channel.

EPICS	EPICS	asyn in-	drvInfo	Description
record	record	terface	string	
name	type			
\$(P)\$(R)	ai	asyn-	TEM-	Temperature input value. This field should be periodically scanned,
		Float64	PERA-	since it is not currently polled in the driver, so I/O Intr scanning
			TURE_IN	_&AhbtHbe used.
\$(P)\$(R)S	cantebbo	asynInt32	TEM-	Temperature scale (units) for this temperature input channel.
			PERA-	Choices are "Celsius" (0), "Fahrenheit" (1), "Kelvin" (2), "Volts"
			TURE_SC	ALE and "Noscale" (5).
\$(P)\$(R)T	Cilijopteo	asynInt32	THER-	Thermocouple type. Choices are "Type J" (1), "Type K" (2), "Type
			MO-	T" (3), "Type 4" (4), "Type R" (5), "Type S" (6), "Type B" (7), "Type
			COU-	N" (8)
			PLE_TYP	E
\$(P)\$(R)F	il <b>ten</b> bbo	asynInt32	TEM-	Temperature filter. Choices are "Filter" (0) and "No filter" (0x400)
			PERA-	
			TURE_FI	LTER

The following is the main medm screen for configuring the analog/temperature inputs on the USB-2408-2AO.

🗙 measCompTemperatureSetup.adl									
	Temp. setup USB2408:								
				Dis	play				
	TC Type	Scale	Filter	Min	Max				
Ti1	Туре К 🗔	Celsius 💷	Filter 💷	20.00	30,00				
Ti2	Туре К 🔲	Celsius 🔳	Filter 💷	0.00	0.00				
Ti3	Туре К 🔲	Celsius 🔳	Filter 💷	0.00	0.00				
Ti4	Туре К 🔲	Fahrenheit 💷	Filter 💷	0.00	0.00				
Ti5	Type J 💷	Celsius 💷	Filter 💷	0.00	0.00				
Ti6	Type J 💷	Celsius 💷	Filter 💷	0.00	0.00				
Ti7	Type J 💷	Celsius 💷	Filter 💷	0.00	0.00				
Ti8	Type J 💷	Celsius 🔲	Filter 💷	0.00	0.00				

Fig. 26: measCompTemperatureSetup.adl

# **Digital I/O Functions**

These are the records defined in the following files:

- measCompBinaryIn.template. This database is loaded once for each binary I/O bit.
- measCompLongIn.template. This database is loaded once for each binary I/O register.
- measCompBinaryOut.template. This database is loaded once for each binary I/O bit.
- measCompLongOut.template. This database is loaded once for each binary I/O register.
- measCompBinaryDir.template. This database is loaded once for each binary I/O bit.

EPICS	EPICS	asyn in-	drvInfo	Description
record	record	terface	string	
name	type			
\$(P)\$(R)	bi	asyn-	DIGI-	Digital input value. The MASK parameter in the INP link defines
		UInt32Dig	ifaAL_INP	UW hich bit is used. The binary inputs are polled by the driver poller
				thread, so these records should have SCAN="I/O Intr".
\$(P)\$(R)	longin	asyn-	DIGI-	Digital input value as a word, rather than individual bits. The MASK
		UInt32Dig	gifaAL_INP	Upparameter in the INP link defines which bits are used. The binary
				inputs are polled by the driver poller thread, so this record should
				have SCAN="I/O Intr".
\$(P)\$(R)	bo	asyn-	DIGI-	Digital output value. The MASK parameter in the INP link defines
		UInt32Dig	itaAL_OU	<b>FRUHT</b> ch bit is used.
\$(P)\$(R)_	RBGiV	asyn-	DIGI-	Digital output value readback. The MASK parameter in the INP link
		UInt32Dig	ifaAL_OU	<b>FRI</b> effnes which bit is used.
\$(P)\$(R)	longout	asyn-	DIGI-	Digital output value as a word, rather than individual bits. The
		UInt32Dig	gittaAL_OU	<b>FRUASK</b> parameter in the INP link defines which bits are used.
\$(P)\$(R)_	RBMgin	asyn-	DIGI-	Digital output value readback as a word, rather than individual bits.
		UInt32Dig	itaAL_OU	<b>PIUSE</b> MASK parameter in the INP link defines which bits are used.
\$(P)\$(R)	bo	asyn-	DIGI-	Direction of this I/O line, "In" (0) or "Out" (1). The MASK param-
		UInt32Dig	ifaAL_DIR	EctelONhe INP link defines which bit is used.

# **Pulse Generator Functions**

Note: These are called "timers" in Measurement Computing's documentation.

These are the records defined in measCompPulseGen.template. This database is loaded once for each pulse generator.

EPICS	EPICS	asyn in-	drvInfo	Description
record	record	terface	string	
name	type			
\$(P)\$(R)R	uho	asyn-	PULSE_R	UNRun" (1) starts the pulse generator, "Stop" (0) stops the pulse gen-
		UInt32		erator. Note that ideally this record should go back to 0 when the
				pulse generator is done, if it is outputting a finite number of pulses
				(see Count record). But unfortunately the Measurement Computing
				library does not have a way to query the status of the timer to see if
				it is done, so this is not possible.
\$(P)\$(R)P	eraiod	asyn-	PULSE_P	ERIGO period, in seconds. The time between pulses can be defined
		Float64		either with the Period or with the Frequency; whenever one record
				is changed the other is updated with the new calculated value.
\$(P)\$(R)F	reaquency	N.A.	N.A.	Pulse frequency, in seconds. The Frequency calculates a new value
				of the Period, and sends the period value to the driver.
\$(P)\$(R)V	Vialoh	asyn-	PULSE_V	Period-
		Float64		15.625 ns).
\$(P)\$(R)E	ekaoy	asyn-	PULSE_D	Hinkiyal pulse delay in seconds after Run is set to 1.
		Float64		
\$(P)\$(R)C	olontgout	asynInt32	PULSE_C	OWINGEr of pulses to output. If the Count is 0 then the pulse generator
				runs continuously until Run is set to 0.
\$(P)\$(R)I	il <b>bS</b> tate	asynInt32	PULSE_I	DIETE Scile Stile of the pulse output line, "Low" (0) or "High" (1). This
				determines the polarity of the pulse, i.e. positive going or negative
				going.

# **Waveform Digitizer Functions**

These records are defined in the following files: - measCompWaveformDig.template. This database is loaded once per module. - measCompWaveformDigN.template. This database is loaded for each digitizer input channel.

EPICS	EPICS	asyn in-	drvInfo	Description	
record	record	terface	string		
$\frac{11}{(P)}$	uhon Provinsta	asynInt32	WAVEDI	<b>WithWerR(1)DWTrfs</b> s to digitize. This cannot be more than the value of	
φ(ε )ψ(ι ()]		u05 mm 52		maxInputPoints that was specified in USB1608GConfig.	
\$(P)\$(R)F	ir <b>st6ba</b> n	asynInt32	WAVEDI	<u><b>BHHRSET</b></u> and <b>HAN</b> digitize. "1" (0) to "8" (7). The database currently	
				assumes differential inputs, so only 8 inputs are available, though	
	1 101	X		this can easily be extended to 16.	
\$(P)\$(R)N	umbhans	asynInt32	WAVEDI	y with the comparis a second state of the second se	
				ferential inputs so only 8 inputs are available though this can easily	
				be extended to 16.	
\$(P)\$(R)T	innæWeF	asyn-	WAVEDI	G_TTIMeBasW waveform. These values are calculated when Dwell or	
	form	Float32Ar	ray	NumPoints are changed. It is typically used as the X-axis in plots.	
\$(P)\$(R)C	u <b>loegtiP</b> oint	asynInt32	WAVEDI	G_COURAFENT performance collected. This does not always increment	
	vvorll	0.0110	WAVEDI	by 1 because the device can transfer data in blocks.	
\$(P)\$(K)L	Watell	asyn- Float64	WAVEDI	<u>j</u> <b>Luevin a</b> per point in seconds. The minimum time is 2 microseconds times NumChans.	
\$(P)\$(R)T	otailTime	asyn-	WAVEDI	G <b>THOTALI</b> ( <b>THMHO</b> digitize NumChans*NumPoints.	
		Float64			
\$(P)\$(R)E	xbFrigger	asynInt32	WAVEDI	G_TEX TriggelG61ERe, "Internal" (0) or "External" (1).	
\$(P)\$(R)E	x <b>t</b> 6lock	asynInt32	WAVEDI	G_TEX TELCCKQGIKce, "Internal" (0) or "External" (1). If External is	
				used then the Dwell record does not control the digitization rate, it is controlled by the external clock. However Dwell should be set	
				to approximately the correct value if possible, because that controls	
				what type of data transfers the device uses.	
\$(P)\$(R)C	o <b>ht</b> inuous	asynInt32	WAVEDI	G_VARNASTENEUCODEs-shot" (0) or "Continuous" (1). This controls	
				whether the device stops when acquisition is complete, or immedi-	
				ately begins another acquisition. Typically "One-shot" is used, be-	
				cause the driver is currently not double-buffered, so data could be overwritten before the driver has a chance to read the data. One ex-	
				ception is when using Retrigger=Enable and TriggerCount less than	
				NumPoints. In that case each trigger will only collect TriggerCount	
				samples, and one wants to use Continuous so that it collects the next	
	1.0.	L (22		TriggerCount samples on the next trigger input.	
(P)(R)	uborRestart	asynInt32	WAVEDI	<b>J_Values</b> <u>Date</u> <b>CS</b> <u>Lisa</u> <u>B</u> <u>B</u> <u>e</u> <sup>*</sup> (U) and "Enable" (1). This controls whether the driver automatically starts another acquire when the pravious one	
				completes. This is different from Continuous mode described above	
				because this is a software restart that only happens after the driver	
				has read the buffer from the previous acquisition.	
\$(P)\$(R)F	etnogger	asynInt32	WAVEDI	G_VARIER KAGHR sable" (0) and "Enable" (1). This controls whether	
			11/11 / PP	the device rearms the trigger input after a trigger is received.	
\$(P)\$(R)T	ri <b>ggeg6ia</b> unt	asynInt32	WAVEDI	J <b>I hist Countrals (DOW) Nit</b> any samples are collected on each trigger input.	
				NumPoints, Retrigger=Enable and Continuous=Enable then each	
				time a trigger is received TriggerCount samples will be collected.	
\$(P)\$(R)E	ubstMode	asynInt32	WAVEDI	G_VAURSafre_MIDIBABle" (0) and "Enable" (1). This controls whether the	
				device digitizes all NumChans channels as quickly as possible dur-	
				ing each sample, or whether it digitizes successive channels at evenly	
				spaced time intervals during the Dwell time. Enabling BurstMode means that all channels are digitized 2 microseconds apart. This can	
				reduce the accuracy if the channels have very different voltages be-	
				cause of the settling time and slew rate limitations of the system.	
\$(P)\$(R)	ubusy	asynInt32	WAVEDI	<b>G_Values</b> are "Stop" (0) and "Run" (1). This starts and stops the wave-	
36	1 11-1			form digitizer.	
<b>~</b> \$(P)\$(R) <b>F</b>	eabodistyF	asynInt32	WAVEDI	j Values Daw FDone" (0) and "Read" (1). "This reads the waveform	
				the driver always reads device when acquisition stops so for quick	
				acquisitions this record can be Passive. To see partial data during	



This is a plot of a digitized waveform captured of someone speaking into a microphone.

Fig. 27: Waveform digitizer plot

#### **Waveform Generator Functions**

These records are defined in the following files: - measCompWaveformGen.template. This database is loaded once per module. - measCompWaveformGenN.template. This database is loaded for each waveform generator output channel.

EPICS	EPICS	asyn in-	drvInfo	Description	
record	record	terface	string		
name	type				
\$(P)\$(R)N	uha Points	asynInt32	WAVE-	Number of points output waveform. The value of this record is equal	
	U	5	GEN NU	MtoPONErINSmPoints if user-defined waveforms are selected, or Int-	
			_	NumPoints if internal predefined waveforms are selected.	
\$(P)\$(R)U	Js <b>toNgonP</b> oi	ntassynInt32	WAVE-	Number of points in user-defined output waveforms. This cannot	
	8		GEN US	ERSENTION PROMINES value of maxOutputPoints that was specified in	
			_	USB1608GConfig.	
\$(P)\$(R)I	nt Normen Retint	s asynInt32	WAVE-	Number of points in internal predefined output waveforms. This can-	
,	U	5	GEN INT	NUL IN TROUNTS the value of maxOutputPoints that was specified in	
			_	USB1608GConfig.	
\$(P)\$(R)U	JsøratieneWI	asyn-	WAVEDI	G Tuse Base IN the form for user-defined waveforms. These values are	
	form	Float32Ar	rav	calculated when UserDwell or UserNumPoints are changed. It is	
				typically used as the X-axis in plots.	
\$(P)\$(R)I	nt Wanne WF	asyn-	WAVE-	Timebase waveform for internal predefined waveforms. These val-	
	form	Float32Ar	ragen int	<b>TEME</b> Walfculated when IntDwell or IntNumPoints are changed. It	
			-	is typically used as the X-axis in plots.	
\$(P)\$(R)C	lu <b>hoenstiP</b> loint	asynInt32	WAVE-	The current point being output. This does not always increment by	
φ(1)φ(1)¢		usj111110 <b>-</b>	GEN CU	<b>RB</b> N <b>CEUROTNCEUCCDCEUCDCDCEUCDCDCEUCDCDCEUCDCDCEUCDCDCEUCDDCDCDCDCDCDCDDCDCDCDCDCDCDCDCDCDCDCDCDCDCDCDCDCDDCDCDDCDCDCDCDDCDCDDCDDCDDCDDDDCDDDDDDDDDDDDD</b>	
\$(P)\$(R)F	regiuency	asvn-	WAVE-	The output frequency (waveforms/second). The value of this record	
φ( <b>i</b> )φ( <b>i</b> ) <b>i</b>	raqueney	Float64	GEN FRI	<b>CON HERVALY</b> UserFrequency if user-defined waveforms are selected or	
		Tioutor		IntFrequency if internal predefined waveforms are selected	
\$(P)\$(R)[	Juaill	asvn-	WAVE-	The output dwell time or period (seconds/sample). The value of this	
ψ(Γ)ψ(Γ()Γ		Float64	GEN DW	Flebord is equal to UserDwell if user-defined waveforms are selected	
		1 100104		or IntDwell if internal predefined waveforms are selected	
\$(P)\$(R)I	[serDwell	asvn-	WAVE-	The output dwell time or period (seconds/sample) for user-defined	
φ(Γ)φ(Γς)ς	Sull well	Float64	GEN US	<b>FreeDwarthals</b> This record is automatically changed if UserFrequency	
		Tioutor		is modified	
\$(P)\$(R)]	ntDwell	asvn-	WAVE-	The output dwell time or period (seconds/sample) for internal pre-	
φ(Γ)φ(Γ)Γ		Eloat64	GEN INT	<b>Dwidd</b> Iwaveforms. This record is automatically changed if IntEre-	
		Tioutor		quency is modified	
\$(P)\$(R)I	Is <b>ør</b> Frequen	CNN A	ΝΑ	The output frequency (waveforms/second) for user-defined wave-	
$\psi(\mathbf{I})\psi(\mathbf{I}\mathbf{V})\mathbf{C}$	sur requen	c yr <b>(</b> .7 <b>r</b> .	1 1.7 1.	forms. This record computes UserDwell and writes to that record	
				This record is automatically changed if UserDwell is modified	
\$(P)\$(R)]	ntForequency	ΝΑ	ΝΔ	The output frequency (waveforms/second) for internal predefined	
φ(Γ)φ(Γ)Γ	litutequency	1 1.7 1.	1 1.7 1.	waveforms. This record computes IntDwell and writes to that record	
				This record is automatically changed if IntDwell is modified	
\$(P)\$(R)7	otailTime	asvn_	WAVE-	The total time to output the waveforms. This is Dwell*NumPoints	
ψ(Γ)ψ(Γ)Γ		Eloat64	GEN TO	TAL TIME	
\$(D)\$(D)F	vbFrigger	asynInt32	WAVE	The trigger source "Internal" (0) or "External" (1)	
φ(1 )Φ(ΙΚ)Ε	- rou igget	asymmu32	GEN EV	T TRIGGER	
\$(D)\$(D)I	vtGlock	asynInt27	WAVE	The clock source "Internal" (1) or "External" (1) If External is used	
φ(Γ)Φ(Κ)Ε	AUDIOCK	asymmit32	GEN EV	The clock source, internal (0) of External (1). If External is used	
			UEN_EA	by the external clock. Howavar Dwall should be set to approximately	
				the correct value if possible, because that controls what tupe of date	
				transfers the device uses	
¢(D)¢(D)¢	ohtinuous	0.01mT-+22	WAVE	Unisters the device uses.	
\$(r)\$(K)(		asymmu32	CEN CO	values are One-shot (0) or Continuous (1). Inis controls	
			GEN_CO	immediately begins again at the start of the survey form	
	atriacer	0.01mT-+20	WAVE	Minimulately begins again at the start of the Waveform.	
\$(P)\$(K)F	enogger	asyn1nt32	WAVE-	values are Disable (U) and Enable (1). This controls whether	
•(D)•(D)	luiber of t	- age 1: (22	UEN_KE	This controls how more enlarge and entry to the second definition of th	
\$(P)\$(K)]	LIRBEROMUN	asynint32	WAVE-	This controls now many values are output on each trigger input.	
			GEN_TR	NumPrints Detainers Failer Fai	
38				time a trigger is received Trigger 2. Table of Contents	
<u> </u>		an 1.00	MAN / F	Unite a ungger is received inggerCount'samples will be output.	
<i>ф</i> (Р)¢(К)ŀ	unusy	asynInt32	WAVE-	values are "Stop" (0) and "Kun" (1). This starts and stops the wave-	
<b>\$</b> ( <b>D</b> ) <b>*</b> ( <b>D</b> ) <b>*</b>			GEN_RU	Niorm generator.	
E \$(P)\$(R)[	IS EN AVEL	asvn-	WAVE-	I use waveform record contains the user-defined waveform generator	



Fig. 28: Plot of an internal predefined waveform (sin wave)

# **Trigger Functions**

These records are defined in measCompTrigger.template. This database is loaded once per module.

EPICS	EPICS	asyn in-	drvInfo	Description
record	record	terface	string	
name	type			
\$(P)\$(R)N	ladbbo	asynInt32	TRIG-	The mode of the external trigger input. Choices are "Positive edge",
			GER_MO	DENegative edge", "High", and "Low".

# 2.2.5 Box for USB-CTR08, USB-3104, and USB-1808X

The following photo is a box we built to house the USB-CTR08, USB-3104, and USB-1808X and provide BNC I/O connections.



Fig. 29: Plot of a user-defined waveform (sum of sin and cos waves)



Fig. 30: GSECARS designed box for USB-CTR08, USB-3104, and USB-1808X

# 2.2.6 Box for USB-2408-2AO

The following photos show a box we built to house the USB-2408-2AO and provide I/O connections.

This is the top view.



Fig. 31: Top view of USB-2408-2AO box

These are the side views.

# 2.2.7 Performance measurements

The following summarizes a simple test of the precision and accuracy of the analog outputs and analog inputs of the USB-1608GX-2AO. The test configuration was with Analog Output 0 connected to Analog Input 0, and also to a Keithley 2700 digital multimeter. The Keithley is a 6.5 digit (22 bit) device, so it can be used to measure the accuracy of the USB-1608GX-2AO analog output, and provide the "true" value to measure the accuracy of the analog input. The 1608GX analog inputs records and the Keithley input had SCAN=0.1 second, so new readings were being made at 10Hz. The following IDL test program was used to drive the analog output from -10V to +10V in 0.1V steps. 10 readings were made of the 1608GX analog inputs, and one reading of the Keithley at each voltage step. These tests were done with the +-10V range of the analog outputs and analog inputs. Since these are 16-bit devices, one bit is 20V/65536 = 0.000305 volts.



Fig. 32: Side views of USB-2408-2AO box

```
pro test_analog_performance_1608, ao=ao, ai=ai, min_volts=min_volts, max_volts=max_
\rightarrow volts, $
                                   step_volts=step_volts, num_samples=num_samples,
→delay=delay, $
                                   keithley=keithley, results
  if (n_elements(ao)
                              eq 0) then ao
                                                     = '1608G:Ao1'
  if (n_elements(ai)
                                                     = '1608G:Ai1'
                              eq 0) then ai
  if (n_elements(min_volts)
                              eq 0) then min_volts = -10.0
  if (n_elements(max_volts)
                              eq 0) then max_volts = 10.0
                              eq 0) then step_volts = 0.1
  if (n_elements(step_volts)
  if (n_elements(num_samples) eq 0) then num_samples = 10
  if (n_elements(delay)
                              eq 0) then delay
                                                     = 0.1
  if (n_elements(keithley)
                              eq 0) then keithley
                                                     = '13LAB:DMM2Dmm_raw.VAL'
  output = min_volts
  samples = dblarr(num_samples)
  num_points = ((max_volts - min_volts) / step_volts + 0.5) + 1
  results = dblarr(4, num_points)
  for i=0, num_points-1 do begin
   output = min_volts + i*step_volts
   t = caput(ao, output)
   wait, 2*delay
   for j=0, num_samples-1 do begin
     wait, delay
     t = caget(ai, temp)
      samples[j] = temp
   endfor
   m = moment(samples)
   results[0,i] = output
   results[1,i] = m[0]
   results[2,i] = sqrt(m[1])
   t = caget(keithley, temp)
   results[3,i] = temp
   print, results[0,i], results[1,i], results[2,i], results[3,i]
  endfor
end
```

The following plot shows the difference of the nominal USB-1608GX-2AO analog output voltage from the Keithley 2700 reading. The mean error is 0.000312V, or just over 1 bit. The RMS error is 0.000203V, or less than 1 bit.

The following plot shows the difference of the mean of 10 readings of the 1608GX analog input voltage from the Keithley 2700 reading. The mean error is 0.000106V, less than 1 bit. The RMS error is 0.000259V, also less than 1 bit.

The following plot shows the standard deviation of 10 readings of the 1608GX analog input voltage. The values range from about 0.001V (~3 bits) at +-10V to less than 0.0003V (1 bit) between -2 and +2V.

The following table contains all of the results from the tests.

1608GX analog output (nominal)	1608GX analog input (mean of 10 readings)	Std. Dev. of 10 1608GX analog input r
-10.00000	-9.99930	0.00084
-9.90000	-9.89978	0.00130
-9.80000	-9.79986	0.00126

1608GX analog output (nominal)	1608GX analog input (mean of 10 readings)	Std. Dev. of 10 1608GX analog input r
-9.70000	-9.69964	0.00134
-9.60000	-9.60018	0.00123
-9.50000	-9.50057	0.00099
-9.40000	-9.40020	0.00117
-9.30000	-9.30010	0.00080
-9.20000	-9.20046	0.00105
-9.10000	-9.09996	0.00118
-9.00000	-9.00035	0.00122
-8.90000	-8.90016	0.00079
-8.80000	-8.80061	0.00118
-8.70000	-8.69996	0.00138
-8.60000	-8.60044	0.00112
-8.50000	-8.50004	0.00098
-8.40000	-8.39973	0.00103
-8.30000	-8.29975	0.00132
-8.20000	-8.19965	0.00108
-8.10000	-8.09986	0.00115
-8.00000	-8.00040	0.00079
-7.90000	-7.90021	0.00088
-7.80000	-7.79950	0.00107
-7.70000	-7.69998	0.00099
-7.60000	-7.60018	0.00092
-7.50000	-7.49990	0.00080
-7.40000	-7.39986	0.00097
-7.30000	-7.29992	0.00101
-7.20000	-7.20006	0.00085
-7.10000	-7.09953	0.00100
-7.00000	-7.00060	0.00088
-6,90000	-6.89986	0.00097
-6.80000	-6.79988	0.00089
-6.70000	-6.69984	0.00107
-6.60000	-6.60017	0.00091
-6.50000	-6.49958	0.00088
-6.40000	-6.40043	0.00105
-6.30000	-6.30005	0.00088
-6.20000	-6.20008	0.00085
-6.10000	-6.10016	0.00076
-6.00000	-6.00052	0.00068
-5.90000	-5.89963	0.00077
-5.80000	-5.80050	0.00076
-5.70000	-5.70013	0.00066
-5.60000	-5.60006	0.00066
-5.50000	-5.50008	0.00082
-5 40000	-5 39989	0.00090
-5.30000	-5.29982	0.00081
-5 20000	-5 19997	0.00087
-5 10000	-5 10021	0.00048
-5.00000	-5 00011	0.00054
-4 90000	-4 89986	0.00071
	110//00	0.000/1

Table 1 – continued from previous page

## measComp

1608GX analog output (nominal)	1608GX analog input (mean of 10 readings)	Std. Dev. of 10 1608GX analog input r
-4.80000	-4.79976	0.00070
-4.70000	-4.69960	0.00082
-4.60000	-4.60090	0.00054
-4.50000	-4.50050	0.00072
-4.40000	-4.40012	0.00076
-4.30000	-4.30039	0.00045
-4.20000	-4.20005	0.00066
-4.10000	-4.10010	0.00068
-4.00000	-4.00012	0.00062
-3.90000	-3.90018	0.00060
-3.80000	-3.80002	0.00059
-3.70000	-3.70019	0.00049
-3.60000	-3.60027	0.00056
-3.50000	-3.50042	0.00063
-3.40000	-3.40017	0.00048
-3.30000	-3.30043	0.00045
-3.20000	-3.20034	0.00064
-3.10000	-3.10027	0.00066
-3.00000	-3.00047	0.00043
-2.90000	-2.90025	0.00060
-2.80000	-2.80021	0.00044
-2.70000	-2.70033	0.00038
-2.60000	-2.60011	0.00058
-2.50000	-2.50001	0.00063
-2.40000	-2.40015	0.00051
-2.30000	-2.29960	0.00043
-2.20000	-2.20050	0.00041
-2.10000	-2.10040	0.00048
-2.00000	-2.00012	0.00054
-1.90000	-1.90018	0.00044
-1.80000	-1.80026	0.00044
-1.70000	-1.70025	0.00062
-1.60000	-1.60043	0.00041
-1.50000	-1.50054	0.00044
-1.40000	-1.40035	0.00037
-1.30000	-1.30001	0.00043
-1.20000	-1.20006	0.00035
-1.10000	-1.10024	0.00048
-1.00000	-1.00035	0.00052
-0.90000	-0.90056	0.00036
-0.80000	-0.80052	0.00050
-0.70000	-0.70011	0.00041
-0.60000	-0.60029	0.00036
-0.50000	-0.50056	0.00035
-0.40000	-0.40031	0.00032
-0.30000	-0.30042	0.00030
-0.20000	-0.20053	0.00048
-0.10000	-0.10037	0.00041
0.00000	0.00018	0.00030

und fr - ntir Tabl 4 .

1608GX analog output (nominal)	1608GX analog input (mean of 10 readings)	Std. Dev. of 10 1608GX analog input re
0.10000	0.09986	0.00046
0.20000	0.19995	0.00032
0.30000	0.30005	0.00035
0.40000	0.39979	0.00046
0.50000	0.49979	0.00032
0.60000	0.60008	0.00028
0.70000	0.69941	0.00041
0.80000	0.79979	0.00019
0.90000	0.89986	0.00037
1.00000	0.99956	0.00032
1.10000	1.09966	0.00051
1.20000	1.19982	0.00045
1.30000	1.29940	0.00041
1.40000	1.39959	0.00041
1.50000	1.49990	0.00035
1.60000	1.59969	0.00035
1.70000	1.69979	0.00052
1.80000	1.80029	0.00016
1.90000	1.89944	0.00050
2.00000	1.99966	0.00047
2.10000	2.09973	0.00045
2.20000	2.19980	0.00041
2.30000	2.29984	0.00044
2.40000	2.40006	0.00023
2.50000	2.49934	0.00032
2.60000	2.59937	0.00038
2.70000	2.69963	0.00054
2.80000	2.79994	0.00032
2.90000	2.90010	0.00033
3.00000	3.00026	0.00021
3.10000	3.09990	0.00027
3.20000	3.19976	0.00041
3.30000	3.30022	0.00022
3.40000	3.39977	0.00061
3,50000	3.49990	0.00045
3.60000	3.59991	0.00068
3.70000	3.69952	0.00039
3.80000	3.79974	0.00052
3.90000	3.89969	0.00043
4,0000	3,99994	0.00029
4 10000	4 09967	0.00042
4 20000	4 19974	0.00063
4.30000	4.29950	0.00058
4 40000	4 39973	0.00066
4 50000	4 50001	0.00055
4 60000	4 60005	0.00048
4 70000	4 70014	0.00043
4 80000	4 70082	0.00059
4 90000	4 80905	0.00059
т.20000	т.07773	0.00002

Table 1 – continued from previous page

## measComp

1608GX analog output (nominal)	1608GX analog input (mean of 10 readings)	Std. Dev. of 10 1608GX analog input r
5.00000	4.99925	0.00059
5.10000	5.09960	0.00066
5.20000	5.19963	0.00087
5.30000	5.29952	0.00072
5.40000	5.39925	0.00084
5.50000	5.49926	0.00059
5.60000	5.59918	0.00065
5.70000	5.70004	0.00073
5.80000	5.79989	0.00081
5.90000	5.89972	0.00087
6.00000	6.00000	0.00076
6.10000	6.10001	0.00038
6.20000	6.19986	0.00047
6.30000	6.29947	0.00071
6.40000	6.39973	0.00077
6.50000	6.49986	0.00068
6.60000	6.60005	0.00091
6.70000	6.69947	0.00085
6.80000	6.79939	0.00065
6.90000	6.89924	0.00083
7.00000	6.99989	0.00074
7.10000	7.09972	0.00091
7.20000	7.20012	0.00074
7.30000	7.30004	0.00073
7.40000	7.39934	0.00061
7.50000	7.50002	0.00073
7.60000	7.60003	0.00074
7.70000	7.69967	0.00101
7.80000	7.79947	0.00089
7.90000	7.89972	0.00094
8.00000	8.00027	0.00083
8.10000	8.09934	0.00090
8.20000	8.19971	0.00095
8.30000	8.29963	0.00112
8.40000	8.39997	0.00073
8.50000	8.49903	0.00089
8.60000	8.59962	0.00080
8.70000	8.69950	0.00109
8.80000	8.79945	0.00084
8.90000	8.89973	0.00111
9.00000	8.99980	0.00083
9.10000	9.09993	0.00071
9.20000	9.19966	0.00098
9.30000	9.29918	0.00090
9.40000	9.39910	0.00097
9.50000	9.49987	0.00106
9.60000	9.59890	0.00102
9.70000	9.70004	0.00110
9.80000	9.79974	0.00105

Table 1 – continued from previous page

1608GX analog output (nominal)	1608GX analog input (mean of 10 readings)	Std. Dev. of 10 1608GX analog input r
9.90000	9.89935	0.00112
10.00000	9.99951	0.00058

Table 1 - continued from previous page

Suggestions and Comments to: Mark Rivers : (rivers@cars.uchicago.edu)

# 2.3 Driver for the USB-CTR08

author

Mark Rivers, University of Chicago

#### Contents

- Driver for the USB-CTR08
  - Introduction
  - Configuration
  - Databases
    - \* Digital I/O Functions
    - \* Pulse Generator Functions
    - \* Scaler Record Support
    - \* Multi-Channel Scaler (MCS) Support
    - \* medm screens
  - Wiring to BCDA BC-020 LEMO Breakout Panels
    - \* Wiring table
  - Performance measurements
  - Restrictions

# 2.3.1 Introduction

This is an EPICS driver for the USB-CTR04 and USB-CTR08 counter/timer modules from MeasurementComputing.

The driver is written in C++, and consists of a class that inherits from asynPortDriver, which is part of the EPICS asyn module.

This module has the following features:

- Digital inputs/outputs
  - 8 signals, individually programmable as inputs or outputs
- Pulse generators. 4 pulse generators each with



Fig. 33: USB-1608GX-2AO analog output voltage error



Fig. 34: USB-1608GX-2AO analog input voltage error



Fig. 35: USB-1608GX-2AO analog input standard deviation



Fig. 36: Photo of USB-CTR08

- 48MHz clock, 32-bit registers
- Programmable period, width, number of pulses, polarity
- Counters. 8 counters (USB-CTR08) or 4 counters (USB-CTR04)
  - 48 MHz maximum count rate
  - Support for EPICS scaler record (similar to Joerger VSC and SIS3820)
  - Support for Multi-Channel Scaler (MCS) mode, similar to SIS3820.

# 2.3.2 Configuration

The following lines are needed in the EPICS startup script for the USBCTR.

```
# This line is for Linux only
cbAddBoard("USB-CTR", "")
## Set the minimum sleep time to 1 ms
asynSetMinTimerPeriod(0.001)
## Configure port driver
# USBCTRConfig(portName,
                               # The name to give to this asyn port driver
               boardNum,
                               # The number of this board assigned by the Measurement
#
→Computing Instacal program
#
               maxTimePoints) # Maximum number of time points for MCS
USBCTRConfig("$(PORT)", 0, 2048, .01)
#asynSetTraceMask($(PORT), 0, TRACE_ERROR|TRACE_FLOW|TRACEIO_DRIVER)
dbLoadTemplate("USBCTR.substitutions")
# This loads the scaler record and supporting records
dbLoadRecords("$(SCALER)/db/scaler.db", "P=USBCTR:, S=scaler1, DTYP=Asyn Scaler.
\rightarrowOUT=@asyn(USBCTR), FREQ=10000000")
# This database provides the support for the MCS functions
dbLoadRecords("$(MEASCOMP)/measCompApp/Db/measCompMCS.template", "P=$(PREFIX), PORT=
\hookrightarrow (PORT)")
# Load either MCA or waveform records below
# The number of records loaded must be the same as MAX_COUNTERS defined above
# Load the MCA records
#dbLoadRecords("$(MCA)/mcaApp/Db/simple_mca.db", "P=$(PREFIX), M=$(RNAME)1, _
→DTYP=asynMCA, INP=@asyn($(PORT) 0), PREC=3, CHANS=$(MAX_POINTS)")
#dbLoadRecords("$(MCA)/mcaApp/Db/simple_mca.db", "P=$(PREFIX), M=$(RNAME)2, _
→DTYP=asynMCA, INP=@asyn($(PORT) 1), PREC=3, CHANS=$(MAX_POINTS)")
#dbLoadRecords("$(MCA)/mcaApp/Db/simple_mca.db", "P=$(PREFIX), M=$(RNAME)3, _
→DTYP=asynMCA, INP=@asyn($(PORT) 2), PREC=3, CHANS=$(MAX_POINTS)")
#dbLoadRecords("$(MCA)/mcaApp/Db/simple_mca.db", "P=$(PREFIX), M=$(RNAME)4, _
→DTYP=asynMCA, INP=@asyn($(PORT) 3), PREC=3, CHANS=$(MAX_POINTS)")
#dbLoadRecords("$(MCA)/mcaApp/Db/simple_mca.db", "P=$(PREFIX), M=$(RNAME)5, _
→DTYP=asynMCA, INP=@asyn($(PORT) 4), PREC=3, CHANS=$(MAX_POINTS)")
```

(continues on next page)

(continued from previous page)

#dbLoadRecords("\$(MCA)/mcaApp/Db/simple\_mca.db", "P=\$(PREFIX), M=\$(RNAME)6, \_ →DTYP=asynMCA, INP=@asyn(\$(PORT) 5), PREC=3, CHANS=\$(MAX\_POINTS)") #dbLoadRecords("\$(MCA)/mcaApp/Db/simple\_mca.db", "P=\$(PREFIX), M=\$(RNAME)7, \_ →DTYP=asynMCA, INP=@asyn(\$(PORT) 6), PREC=3, CHANS=\$(MAX\_POINTS)") #dbLoadRecords("\$(MCA)/mcaApp/Db/simple\_mca.db", "P=\$(PREFIX), M=\$(RNAME)8, \_ →DTYP=asynMCA, INP=@asyn(\$(PORT) 7), PREC=3, CHANS=\$(MAX\_POINTS)") #dbLoadRecords("\$(MCA)/mcaApp/Db/simple\_mca.db", "P=\$(PREFIX), M=\$(RNAME)9, \_ →DTYP=asynMCA, INP=@asyn(\$(PORT) 8), PREC=3, CHANS=\$(MAX\_POINTS)") # This loads the waveform records dbLoadRecords("\$(MCA)/mcaApp/Db/SIS38XX\_waveform.template", "P=\$(PREFIX), R=\$(RNAME)1,  $\rightarrow$  INP=@asyn(\$(PORT) 0), CHANS=\$(MAX\_POINTS)") dbLoadRecords("\$(MCA)/mcaApp/Db/SIS38XX\_waveform.template", "P=\$(PREFIX), R=\$(RNAME)2, →INP=@asyn(\$(PORT) 1), CHANS=\$(MAX\_POINTS)") dbLoadRecords("\$(MCA)/mcaApp/Db/SIS38XX\_waveform.template", "P=\$(PREFIX), R=\$(RNAME)3, →INP=@asyn(\$(PORT) 2), CHANS=\$(MAX\_POINTS)") dbLoadRecords("\$(MCA)/mcaApp/Db/SIS38XX\_waveform.template", "P=\$(PREFIX), R=\$(RNAME)4, →INP=@asyn(\$(PORT) 3), CHANS=\$(MAX\_POINTS)") dbLoadRecords("\$(MCA)/mcaApp/Db/SIS38XX\_waveform.template", "P=\$(PREFIX), R=\$(RNAME)5, \_ →INP=@asyn(\$(PORT) 4), CHANS=\$(MAX\_POINTS)") dbLoadRecords("\$(MCA)/mcaApp/Db/SIS38XX\_waveform.template", "P=\$(PREFIX), R=\$(RNAME)6, →INP=@asyn(\$(PORT) 5), CHANS=\$(MAX\_POINTS)") dbLoadRecords("\$(MCA)/mcaApp/Db/SIS38XX\_waveform.template", "P=\$(PREFIX), R=\$(RNAME)7, →INP=@asyn(\$(PORT) 6), CHANS=\$(MAX\_POINTS)") dbLoadRecords("\$(MCA)/mcaApp/Db/SIS38XX\_waveform.template", "P=\$(PREFIX), R=\$(RNAME)8, →INP=@asyn(\$(PORT) 7), CHANS=\$(MAX\_POINTS)") dbLoadRecords("\$(MCA)/mcaApp/Db/SIS38XX\_waveform.template", "P=\$(PREFIX), R=\$(RNAME)9, →INP=@asyn(\$(PORT) 8), CHANS=\$(MAX\_POINTS)") asynSetTraceIOMask(\$(PORT),0,2) #asynSetTraceFile("\$(PORT)",0,"\$(MODEL).out") < save restore.cmd save\_restoreSet\_status\_prefix(\$(PREFIX)) dbLoadRecords("\$(AUTOSAVE)/asApp/Db/save\_restoreStatus.db", "P=\$(PREFIX)") iocInit seq(USBCTR\_SNL, "P=\$(PREFIX), R=\$(RNAME), NUM\_COUNTERS=\$(MAX\_COUNTERS), FIELD=\$(FIELD)") create\_monitor\_set("auto\_settings.req",30)

The measComp module comes with an example iocBoot/iocUSBCTR directory that contains and example startup script and example substitution files.

# 2.3.3 Databases

The following tables list the database template files that are used with the USB-CTR04/08.

# **Digital I/O Functions**

These are the records defined in the following files:

- measCompBinaryIn.template. This database is loaded once for each binary I/O bit.
- measCompLongIn.template. This database is loaded once for each binary I/O register.
- measCompBinaryOut.template. This database is loaded once for each binary I/O bit.
- measCompLongOut.template. This database is loaded once for each binary I/O register.
- measCompBinaryDir.template. This database is loaded once for each binary I/O bit.

EPICS	EPICS	asyn in-	drvInfo	Description
record	record	terface	string	
name	type		_	
\$(P)\$(R)	bi	asyn-	DIGI-	Digital input value. The MASK parameter in the INP link defines
		UInt32Dig	gital_INP	UW hich bit is used. The binary inputs are polled by the driver poller
				thread, so these records should have SCAN="I/O Intr".
\$(P)\$(R)	longin	asyn-	DIGI-	Digital input value as a word, rather than individual bits. The MASK
		UInt32Dig	itaaL_INP	Upparameter in the INP link defines which bits are used. The binary
				inputs are polled by the driver poller thread, so this record should
				have SCAN="I/O Intr".
\$(P)\$(R)	bo	asyn-	DIGI-	Digital output value. The MASK parameter in the INP link defines
		UInt32Dig	itaaL_OU	<b>FRUM</b> Tch bit is used.
\$(P)\$(R)_	RBGiV	asyn-	DIGI-	Digital output value readback. The MASK parameter in the INP link
		UInt32Dig	ifaal_OU	<b>FRI</b> effnes which bit is used.
\$(P)\$(R)	longout	asyn-	DIGI-	Digital output value as a word, rather than individual bits. The
		UInt32Dig	itaaL_OU	<b>FRUASK</b> parameter in the INP link defines which bits are used.
\$(P)\$(R)_	RBMgin	asyn-	DIGI-	Digital output value readback as a word, rather than individual bits.
		UInt32Dig	itaaL_OU	<b>FPUGE</b> MASK parameter in the INP link defines which bits are used.
\$(P)\$(R)	bo	asyn-	DIGI-	Direction of this I/O line, "In" (0) or "Out" (1). The MASK param-
		UInt32Dig	itaaL_DIR	EctalONhe INP link defines which bit is used.

#### **Pulse Generator Functions**

Note: These are called "timers" in Measurement Computing's documentation.

These are the records defined in measCompPulseGen.template. This database is loaded once for each pulse generator.

EPICS	EPICS	asyn in-	drvInfo	Description
record	record	terface	string	
name	type		Ũ	
\$(P)\$(R)F	uho	asyn-	PULSE_R	UNRun" (1) starts the pulse generator, "Stop" (0) stops the pulse gen-
		UInt32		erator. Note that ideally this record should go back to 0 when the
				pulse generator is done, if it is outputting a finite number of pulses
				(see Count record). But unfortunately the Measurement Computing
				library does not have a way to query the status of the timer to see if
				it is done, so this is not possible.
\$(P)\$(R)F	eraiod	asyn-	PULSE_P	ERIGD period, in seconds. The time between pulses can be defined
		Float64		either with the Period or with the Frequency; whenever one record
				is changed the other is updated with the new calculated value.
\$(P)\$(R)F	reaquency	N.A.	N.A.	Pulse frequency, in seconds. The Frequency calculates a new value
				of the Period, and sends the period value to the driver.
\$(P)\$(R)V	Vialoh	asyn-	PULSE_V	Ported width, in seconds. The allowed range is 15.625 ns to (Period-
		Float64		15.625 ns).
\$(P)\$(R)I	)elaay	asyn-	PULSE_D	Elhaval pulse delay in seconds after Run is set to 1.
		Float64		
\$(P)\$(R)C	olontgout	asynInt32	PULSE_C	CNUMBER of pulses to output. If the Count is 0 then the pulse generator
				runs continuously until Run is set to 0.
\$(P)\$(R)I	dl <b>bS</b> tate	asynInt32	PULSE_I	DILFIE Scille Stille Stille the pulse output line, "Low" (0) or "High" (1). This
				determines the polarity of the pulse, i.e. positive going or negative
				going.

# **Scaler Record Support**

The USBCTR driver supports the EPICS scaler record via the devScalerAsyn.c device support originally from the synApps std module but which has been moved into the scaler module. It supports up to 8 channels. The following wiring connections must be made in order for counters 1-8 to be stopped by counter 0, as is normally desired.

• Counter 0 Output must be connected to the Gate input on Counters 1-7.

The .PR1 preset is performed in hardware via the Counter 0 Output and Counters 1-7 gates. Counters 1-7 can also be set as preset counters, and the scaler record will stop counting when any of these preset values (.PR2-.PR8) are exceeded. However, unlike the .PR1 preset, these presets are done in software in the driver polling routine. The device sends readings at 100 Hz, and whenever a preset is exceeded counting is stopped. Each of the counters will have counted for exactly the same amount of time, but the actual count time could be up to 0.01 seconds longer than the time when the preset was reached.

Counter 0 is normally used as the preset counter, and is connected to a fixed frequency source. Any of the on-board pulse generators can be used to provide this frequency source, for example. It is important to set the scaler record .FREQ field to be the value of the Frequency\_RBV of the pulse generator (the actual frequency) and not the Frequency field (the requested frequency) since these can differ, particularly at frequencies >1 MHz.

## Multi-Channel Scaler (MCS) Support

The USBCTR driver provides multi-channel scaler support very similar to the SIS3820 driver in the synApps mca module. The support has the following properties:

- The number of counters being used in MCS mode can be selected with the FirstCounter and LastCounter records. Each can range from 0 to 7; LastCounter must be greater than or equal to FirstCounter. The number of active counters can thus range from 1 to 8.
- The minimum dwell time, either with internal or external channel advance, is 250 ns times the number of active counters. For example if only 2 counters are being used, the clock input on Counter 0 and a signal on Counter 1, then the minimum dwell time is 500 ns. If all 8 counters are being used then the minimum dwell time is 2 microseconds.
- Either MCS or waveform records can be used to hold the time series data.
- There is no limitation on the length of the waveform or mca records, only the size of system RAM.
- An external channel advance signal can be used directly by connecting it to the External Clock Input (CLKI) on the USB-CTR module. The minimum dwell time (period) of this signal is described above.
- An external channel advance can be "prescaled" (frequency divided by N) by connecting it to a counter input. This counter is assigned to the PrescaleCounter record. The Counter Output of the PrescaleCounter must be connected to the External Clock Input on the USB-CTR module. I have asked Measurment Computing to consider adding a prescale register for the CLKI signal in a future firmware version, but I don't know if this will be done.
- To achieve the shortest dwell times the counter must be read in 16-bit mode rather than 32-bit mode. This is handled automatically by the driver. If the dwell time is less than 100 microseconds the counters are read in 16-bit mode, while for longer dwell times they are read in 32-bit mode. There is no possible loss of data when reading in 16-bit mode because at the maximum count rate of 48 MHz only 4800 counts can occur in 100 microseconds, which is much less than the 16-bit limit. NOTE: When using external channel advance the Dwell record should be set to the approximate time between external pulses. This will cause the correct 32-bit/16-bit switch to occur so that the minimum dwell time can be reached and so the counters don't overflow 16-bits for longer dwell times.

The following record are defined in measCompMCS.template. This database is loaded once per module.

EPICS	EPICS record	asyn in-	drvInfo string	Description	
name	type	lenace	Sung		
\$(P)\$(R)S	Nbi_Connec	cted.A.	N.A.	This record is 1 ("Connected") if all PVs have connected in the US-BCTR_SNL State Notation Language program.	
(P)(R)E	rabseAll	asynInt32	MCA_ER	ABEases the MCS data, setting the arrays and the elapsed times to 0.	
\$(P)\$(R)E	rabseStart	asynInt32	MCA_ER	ABFases the MCS data and then starts MCS acquisition by forward linking to StartAll.	
\$(P)\$(R)S	tabtAll	asynInt32	MCA_ST	A BTarts COULS REquisition.	
\$(P)\$(R)A	c <b>quiry</b> ing	N.A.	N.A.	Busy record is 1 ("Acquiring") when MCS is acquiring and 0 ("Done") when done	
\$(P)\$(R)S	tchpoAll	asynInt32	MCA_ST	OBtops QUAR Ecquisition.	
\$(P)\$(R)P	ressortReal	asyn- Float64	MCA_PR	E <b>\$FES</b> eRteAlItime. If non-zero acquisition will stop after this time.	
\$(P)\$(R)E	lapisedReal	asyn- Float64	MCA_EL	AESEBECREALtime.	
\$(P)\$(R)R	ebodAll	N.A	N.A.	Forces a read of all of the array data. This is done by the SNL pro- gram.	
\$(P)\$(R)N	lukonAglobut	asynInt32	MCA_NU	MThehan Nove bestime points to acquire.	
\$(P)\$(R)C	ulocagichan	nælsynInt32	MCS_CU	RIFENELT ROUNTHE point in the acquisition.	
\$(P)\$(R)E	Waell	asyn- Float64	MCA_DW	ETHe_dtvMIR ime per time point in internal channel advance mode.	
\$(P)\$(R)C	h <b>bo</b> nelAdva	anæsynInt32	MCA_CH	<b>ADV</b> <u>c</u> <b>StobBRGE</b> vance source. 0="Internal" uses DWELL record, 1="External" uses External Clock Input on USB-CTR module.	
\$(P)\$(R)P	retssoale	asynInt32	MCA_PR	<b>ESUGAple</b> scale factor for the external channel advance source. To use Prescale the external clock must be input to the counter channel se- lected by PrescaleCounter, and the output of the PrescaleCounter counter channel must be connected to the External Clock Input. Note that due to hardware limitations Prescale must be & gt; 1. For no prescaling the external channel advance source must be connected	
\$( <b>D</b> )\$( <b>D</b> ))	MESCounte	rNGmIht27	ΝΛ	directly to the External Clock Input.	
(N=1-8)	viewscounte	1 las yillatti O Z	<b>N.A.</b>	(1).	
\$(P)\$(R))]	MGCSDIOEn	a <b>ble</b> ynInt32	N.A.	Enable collecting digital I/O word in MCS mode. Choices are "No" (0) and "Yes" (1).	
\$(P)\$(R)P	ra <b>sıbbe</b> Cour	nt <b>es</b> ynInt32	MCS_PRI	ESCHALCE <u>u</u> <b>ATER</b> to use for prescaling the external channel advance in MCS mode. 0="CNTR0" 7="CNTR7".	
\$(P)\$(R)P	oimul6bAcction	asynInt32	MCS_PO	<ul> <li>NI@htAISTION the first time point in the MCS scan is handled. The USB-CTR always reads the current scaler counts as soon as MCS acquisition begins, rather than after the first channel advance occurs. This record selects one of the following 3 modes:</li> <li>"Clear" (0) In this mode the scalers are cleared to 0 before they are read. This means that the counts in first time point for each counter will be 0.</li> <li>"No clear" (1) In this mode the scalers are not cleared before they are read. This means that there will normally be a large number of counts in the first time point, since the counters will have been counting since they were last cleared.</li> <li>"Skip" (2) In this mode the first time point will be skipped, i.e. not read into the mca or waveform records. The first time point will thus contain the counts after MCS acquisition was started until the first channel advance signal is received, either internal or external. This is probably the mode that will be most useful. However, it does require N+1 channel advance signals</li> </ul>	
58				rather than IN. This is handled by the driver for internal channel advance. But for external channel advance the user must en- sure that N+1 pulses are sent. For example if NUseAll=2000 then 2001 pulses must be sent before acquisition will stop.	

#### medm screens

The following is the main medm screen for controlling the USB-CTR04/08.

The following is the medm screen for the EPICS scaler record using the USB-CTR04/08.

The following is the medm screen for controlling the MCS mode of the USB-CTR04/08.

# 2.3.4 Wiring to BCDA BC-020 LEMO Breakout Panels

The following photos show the BCDA BC-020 LEMO breakout panels wired to the USB-CTR08. A BC-020 with a BC-087 daughter card (left) is used for the 8 counter signals, and a BC-020 with wire-wrapping (right) is used for digital I/O, timer output, clock I/O, etc. .

#### Wiring table

Digital	I/O <b>and</b> other signa	als using wire-	wrap connections
50-pin ribbon	USB-1608GX	BC-020	EPICS Function
connector pin	screw terminal	connector	
1	DIOO	J1	Digital I∕O bit 0
2	GND	J1 shell	Ground
3	DI01	J2	Digital I/O bit 1
4	GND	J2 shell	Ground
5	DIO2	J3	Digital I/O bit 2
6	GND	J3 shell	Ground
7	DIO3	]4	Digital I/O bit 3
8	GND	J4 shell	Ground
9	DIO4	J 5	Digital I/O bit 4
10	GND	J5 shell	Ground
11	DIO5	J6	Digital I/O bit 5
12	GND	J6 shell	Ground
13	DIO6	J7	Digital I/O bit 6
14	GND	J7 shell	Ground
15	DIO7	J8	Digital I/O bit 7
16	GND	J8 shell	Ground
17	TMR0	39	Pulse generator 🛛 output
18	GND	J9 shell	Ground
19	TMR1	J10	Pulse generator 1 output
20	GND	J10 shell	Ground
21	TMR2	J11	Pulse generator 2 output
22	GND	J11 shell	Ground
23	TMR3	J12	Pulse generator 3 output
24	GND	J12 shell	Ground
25	TRIG	J13	Trigger input <b>for</b> MCS
26	GND	J13 shell	Ground
27	CLKI	J14	External channel advance input
28	GND	J14 shell	Ground
29	CLK0	J15	Clock output
30	GND	J15 shell	Ground
31	+V0	J16	+5 volt output
32	GND	J16 shell	Ground

(continues on next page)

X USBCTR.adl@corvette			- 🗆 X
USB-CTR08	3 13IDC:U	JSBCTR:	
Model USB-CTR08 UL version 1.2.0			
Model # 295	Driver	version 4.2	2
Unique ID 0214D588	Poll sleep	time (ms)	50.0
Firmware 0.10	Poll cycle	time (ms)	50.1
Dig	ital I/O		
1 2 3	4 5 6	78	
Inputs 🛑 🕘 !	• • •	🕘 🛑 🛛	:11
Outputs Low Low Low I	.ow Low Low	Low Low Dx	0
In In In	In In In	In In	
Direction Out Out Out (	Out Out Out	Out Out	
Pulse	e generato	ors	
9,6000e+06	1,0000e+06	1,0000e+05	100.0000
Frequency 9.6000e+06	1.0000e+06	1.0000e+05	100.0000
1.041/e=0/	1,0000e-06	1,0000e-05	1.0000e-02
Period 1.041/e-0/	1.0000e-06	1.0000e-05	1.0000e-02
0,5000	0,000	0.000	0,000
Duty cycle p.5000	p.5000	p.5000	p.5000
0, 20000-00	5,0000-07	5,0000-06	5,0000-03
0.0000e+00	0.00000000	0.0000e+00	0.0000e+00
Initial delau 0.0000e+00	0.0000e+00	0.0000e+00	0.0000e+00
# pulses	0	0	0
Idle state Low =	Low 🖬	Low 🖃	Low I
Running	Running	Running	Running
Start	Start	Start	Start
Stop	Stop	Stop	Stop
Scaler	MCS	Asvn	record
DScaler	THOSE 1		Basyn

Fig. 37: USBCTR.adl

💐 scaler_full.adl@corvette	_	
Done OneShot time Count time	Elapsed time	Counts
Count AutoCount 2.00 10.000	3.030	Cts/sec
# Description Gate? Preset count	Actual count	Calc result
1. NY 1000000	3030000	1000000,000
2. N Y 10000	303000	100000,000
3. N Y 0	0	0.000
4. <u>NY</u> 0	0	0.000
5. NY 0	0	0.000
6. NY 0	0	0.000
7. NY 0	0	0.000
8. NY 0	3030000	999999, 000
Delay 0.000 (s) Clock 1.000	D <mark>e+06</mark> Hz Displa	yFreq 10.00 Hz
AutoCount: Delay <mark>0.500 (</mark> s)	Displa	yFreq 10,00Hz
Calculations ENABLE =		Less More

Fig. 38: scaler\_full.adl



# Fig. 39: USBCTR\_MCS.adl



Fig. 40: USBCTR\_MCS\_plots.adl



Fig. 41: BC-020 LEMO breakout panels with USBCTR-08



Fig. 42: Top view of USBCTR-08 with BC-020 LEMO breakout panels

(continued from previous page)

Counte	er I/O using wire-w	wrap connection	ns
50-pin ribbon	USB-CTR08	BC-020 EPT	CS Function
connector pin	screw terminal	connector	
1	COTN	11	Scaler 1 input
2	GND	]1 shell	Ground
3	COGT	12	Scaler 1 gate input
4	GND	12 shell	Ground
5	COO	13	Scaler 1 output
6	GND	13 shell	Ground
7	CITN	14	Scaler 2 input
8	GND	J4 shell	Ground
9	C1GT	35	Scaler 2 gate input
10	GND	J5 shell	Ground
11	C10	J6	Scaler 2 output
12	GND	J6 shell	Ground
13	C2IN	J7	Scaler 3 input
14	GND	J7 shell	Ground
15	C2GT	J8	Scaler 3 gate input
16	GND	J8 shell	Ground
17	C20	]9	Scaler 3 output
18	GND	J9 shell	Ground
19	C3IN	J10	Scaler 4 input
20	GND	J10 shell	Ground
21	C3GT	J11	Scaler 4 gate input
22	GND	J11 shell	Ground
23	C40	J12	Scaler 4 output
24	GND	J12 shell	Ground
25	C4IN	J13	Scaler 5 input
26	GND	J14 shell	Ground
27	C4GT	J14	Scaler 5 gate input
28	GND	J14 shell	Ground
29	C40	J15	Scaler 5 output
30	GND	J15 shell	Ground
31	C5IN	J16	Scaler 6 input
32	GND	J16 shell	Ground
33	C5GT	J17	Scaler 6 gate input
34	GND	J17 shell	Ground
35	C50	J18	Scaler 6 output
36	GND	J18 shell	Ground
37	C6IN	J19	Scaler 7 input
38	GND	J19 shell	Ground
39	C6GT	J20 Sca	aler 7 gate input
40	GND	J20 shell	Ground
41	C60	J21	Scaler 7 output
42	GND	J21 shell	Ground
43	C7IN	J22	Scaler 8 input
44	GND	J22 shell	Ground
45	C7GT	J23	Scaler 8 gate input
46	GND	J23 shell	Ground

(continues on next page)

(continued from previous page)

47	C70	J24	Scaler 8 output
48	GND	J24 shell	Ground
In addition to the inputs of counter This is choosen	nese connections c rs 1-7 (C1GT - C7G	ounter 0 outp T) at the mode	ut (C00) was connected to the gate ule screw terminals.

# 2.3.5 Performance measurements

The binary input bits are polled at 100 Hz, and the input records have SCAN=I/O Intr. There is thus a worse-case latency of 0.01 seconds in detecting a transition on these bits.

If the scaler record is run under the following conditions:

- Counter 0 Output connected to the Gate Input of Counters 1-7
- Pulse generator 0 frequency=32 MHz, connected to Counter 0 input
- Pulse generator 1 frequency=32 MHz, connected to Counter 1 input
- Pulse generator 2 frequency=32 MHz, connected to Counter 2 input
- Pulse generator 3 frequency=32 MHz, connected to Counter 3 input
- Scaler record .FREQ field = 3.2e7
- Scaler record preset time = 1.0 second
- Only scaler channel 1 is preset (.G1=Y, .G2-.G8=N)

After each count cycle .S1=32000000 counts exactly, .S2-.S4=32000000 += 1 count. There is thus no cross-talk with all channels running at 32 MHz, and the gate signals are working as designed.

If Pulse Generator 2 is changed to 3.2 MHz, .PR2 is set to 1600000, and .G2 is set to Y, then the scaler is stopped by channel 2 in the software polling routine. In this case it counts for exactly 0.50 seconds. However, if .PR2 is increased to 1600001 then it counts for 0.51 seconds. This corresponds to the worst case error due to the 100 Hz rate at which the scaler values are read. Note that all counters are active for exactly 0.51 seconds, so the counts all accurately reflect this count time. The count time is just slightly longer than requested due to the finite polling interval.

In MCS mode the measured minimum dwell time in both internal and external channel advance mode agrees with the datasheet, i.e. 250 ns \* number of active counters. I was not able to measure any dead time between time bins in MCS mode. When sending exactly 8000000 pulses at 8 MHz to channel 0 with a 1 ms internal dwell time the total number of counts in the MCA record was 8000000. This means that no pulses were lost during the 1000 channel advances that happened during this time.

# 2.3.6 Restrictions

- The EPICS driver only uses the Totalize mode of the counters. With the scaler record it does a one-shot totalize, while in the MCS mode it totalizes into time-bins. The USB-CTR08 is also capable of running in 3 other modes.
  - 1. In Period mode it measures the time between the rising or falling edges of successive input pulses.
  - 2. In Pulse Width measurement mode it measures the time between the rising and falling edges of a each pulse.
  - In Timing Mode it measures the time between an event on the counter input and another event on the counter gate.

None of these modes are currently supported by the EPICS driver, but they could be added in a future release.

- In Totalize mode each counter has many options in how it works: count up/down, gate clears counter, gate controls counter direction, preset counts where the output signal goes high/low, polarity of the output, etc. These options are not currently exposed in the EPICS driver.
- The EPICS driver only works in 32-bit counter depth mode. The USB-CTR08 can count with a 64-bit counter depth. asyn does not currently have support for 64-bit integer data types, so this cannot be supported.
- To work with the scaler record the counter 0 output must be wired to the gate inputs of counters 1-7 as discussed above.

Suggestions and Comments to: Mark Rivers : (rivers@cars.uchicago.edu)